

Baltimore, Maryland



Noncoronary MDCT – Topics:

LV mass

Pericardial disease Congenital disease Pulmonary vein anatomy LV function Valve analysis Myocardial infarction, perfusion







Temporal Resolution Determined by Gantry Rotation Time for LV cines

- Rotation speed: 330-500msec
 - *1/2 scan:* 165-250 msec temporal resolution









- MRI:
- Echo, Nuclear: 25 msec or less



LV function

- Sophisticated segmentation algorithms are adaptive to heart rate, chosen 'automatically' (little user control) to optimize CT angiography.
- Higher segmentation factor beneficial for LV function evaluation.
- Raw data reconstruction available on some scanners; higher segmentation may improve temporal resolution.



LV function: slice thickness

- Instead of 0.5 mm slices for CTA, use 5 mm slices (1000 images).
- MRI uses 6-8 mm slices at 1 cm intervals, ~ 250 images.



MRI: 8 mm slice thickness



LV function – reformat

- MDCT images are acquired in the axial plane
- LV quantitative analysis usually performed in the short axis plane.



Axial images







Papillary Muscles • Correlated with LV wall mass

• Papillary muscle mass accounts for 8.9% of the total LV mass in both men and women

(r=0.81, p<0.001)*

J. Vogel-Claussen, JHU









C	ata courtesy of Fujita Health Uni	versity, Aic	hi, Japan
	LEFT VENTRICULA	R VOLU	VE
	RESULTS		
	Body Surface Area:	1.89	m²
On IOn IOn IOn	ED volume:	357.65	ml
	ED volume/BSA:	189.04	ml/m²
	ES volume:	241.32	ml
AF AF AF AF	ES volume/BSA:	127.55	ml/m²
les her her her	Stroke volume:	116.33	ml
	Stroke volume/BSA:	61.49	ml/m²
C 1400 C 1400 C 1400 C 1400	Ejection fraction:	32.53	%
	LV mass ED:	175.24	g
	LV mass ED/BSA:	92.62	g/m²
	LV mass ES:	190.67	g
A AND A AND A AND AND AND	LV mass ES/BSA:	100.78	g/m²
and the and the and	PER:	281.68	ml/s
	PER/EDV:	0.79	EDV/s
	TPER:	400.00	ms
	TPER phase number:	5	
	PFR:	203.55	ml/s
	PFR/EDV:	0.57	EDV/s
	TPFR:	300.00	ms
	TPFR phase number:	11	
	<u>L</u>		

















Kal Uwe Juergens, MD Matthias Grude, MD David Maintz, MD Eva Maria Falenberg, MD Thomas Wichter, MD, FESC Walter Heindel, MD Roman Fischbach, MD

Index to

Data courtesy of Fujita Health University, Aichi, Japan

Multi-Detector Row CT of Left Ventricular Function with Dedicated Analysis Software versus MR Imaging: Initial Experience¹

- 4 slice scanner
- Temporal resolution 125-250 msec

• Comparable volumes, function compared to 32 msec temporal resolution MRI





Parameter	Row CT*	Imaging*	r Value	P Value
nd-diastolic volume (mL)	104.8 ± 18.7	106.2 ± 19.6	0.98	>.05
nd-systolic volume (mL)	57.2 ± 14.7	57.1 ± 14.5	0.99	>.05
Stroke volume (mL)	47.6 ± 7.4	49.1 ± 7.9	0.92	>.05
jection fraction (%)	46.1 ± 6.5	46.8 ± 5.9	0.97	>.05
Avocardial mass (g)	108 7 + 23 7	109.6 + 27.1	0.92	> 05
PER (mL/sec)	292.0 ± 124.5	229.4 ± 37.3	0.46	.04
Firme to RER (meac)	408.1 ± 180.5	331.0 ± 109.2 122.9 ± 49.2	0.67	05
Time from end systole to PER (msec)	282.9 ± 204.7	155.5 ± 83.4	0.27	.03
† Values are the mean \pm standard dev † Values were determined with the †	iation. t test.	inforat model bo	ort roto 72	106



Why Evaluate Valve Function on MDCT?

Echo: AVA is not determined directly but calculated using the continuity equation.

• TEE and TTE both operator dependent.

• Underestimation of severity due to failure to obtain a parallel intercept angle between the Doppler beam and aortic jet.

• MDCT: Noninvasive, nonoperator-dependent technique for direct measurement of AVA.











MDCT: Assessment of Mitral Valve



Compared to Echo: • Excellent correlation with valve leaflet thickness

• Excellent agreement with Mitral Annulus Calcification

> Figure 3 Willman, J, et al. Eur Radiology 2002; 12

IN	N
lve leaflets	
4 (25)	4 (25)
12 (75)	12 (75)
leification	
7 (78)	6 (67)
2 (22)	3 (33)
	lve leaflets 4 (25) 12 (75) lcification 7 (78) 2 (22)









MDCT Perfusion and Viability Imaging

MDCT Perfusion/ Viability Imaging Rationale:

- Potential to perform CTA, function and viability in a single 15-20 minute exam
- AICD, pacemakers, MRI contraindication
- High spatial resolution (0.4 mm) compared to 6-8 mm slice resolution for MRI.

adapted from A. Lardo

<u>Acute</u> Myocardial Infarction: Contrastenhanced MDCT in a Porcine Model

Udo Hoffman, Ryan Millea, Christian Enzweiler, Maros Ferencik, Scott Gulick, Jim Titus, Stephan Achenbach, Dylan Kwait, David Sosnovik, Thomas J. Brady Radiology, 231:697-701, 2004.

- Porcine AMI model (N=7)
- 4 slice MDCT scanning
- 3 hours post-coronary ligation LAD or LADD
- CT Infarct size 17 ± 6 % similar to TTC 14 ± 6 %



Characterization of <u>Acute</u> MI Using Contrast Enhanced MDCT

Albert C. Lardo, Marco Cordeiro, Veronica Fernandes, Andre Schmidt, Menekhem Zviman, Joao A. C. Lima

Circulation 2004;110 (Supplement):III-522

- Canine AMI model (N=7)
- 32 slice MDCT 1 hour post-reperfusion.
- Imaging each 5 min up to 40 min post-contrast
- CT infarct size 24 ± 7 %, mean difference 4%
- r = 0.93 with TTC

adapted from A. Lardo



















Cardiac Ablation for Arrhythmia

- Catheter based treatment to kill electrically active viable cell islets
- Pre-procedural MDCT for anatomic correlation
- Visualization of RF ablation lesions





MDCT Lesion Morphology



• Hypoenhanced core: coagulative necrosis and microvascular obstruction (contrast does not enter local microcirculation)

Hyperintense periphery: edema, cell necrosis

adapted from A. Lardo

MDCT Perfusion Imaging

Requirements:

- Vasodilator (adenosine)
- During first-pass, contrast-enhanced MDCT.
- Rapid imaging from base to apex.

Possible Advantages:

- Simultaneous coronary imaging.
- Imaging of the entire LV

adapted from A. Lardo



Helical MDCT perfusion imaging in a canine model of LAD stenosis during adenosine infusion. (Gantry rotation time: 400 ms, Detector collimation: 0.5 X 32, tube current: 400 mA, tube voltage: 120 kV, Visipaque[™] 2.5 ml/sec for 100 ml)

George RT, et al. AHA, 2005. adapted from A. Lardo

MDCT Dynamic Perfusion



Serial imaging of the mid-left ventricle over time in a canine model of LAD stenosis. (Detector collimation: 8 mm X 4, 120 kV, 150 mA, VisipaqueTM 10 ml/sec for 3 seconds.)

George RT, et al. Work in progres adapted from A. Larc



Conclusions: Noncoronary MDCT

- Global and regional left ventricular function can be assessed by MDCT and coupled with coronary CTA.
- MDCT stress perfusion, viability and scar imaging has high spatial resolution, promising tool.

Acknowledgements

- Al Lardo, Ph.D.
- Rich George, M.D.
- Jun Dong, M.D.
- Elliot Fishman, M.D.
- João Lima, MD
- Rob van der Geest, PhD

Conclusion - MDCT

Potential for *comprehensive* morphologic and functional assessment.



